## Gebrauchsanweisung Dental Blank

Material:

## Polyan Plus® Economic Line

### 1. Kontaktdaten

 DentalPlus GmbH  
Kohlgrub 5  
83122 Samerberg  
Tel.: 08032 9892007  
Fax: 08032 9882790  
E- Mail: info@dentalplus.info  
Internet: www.dentalplus.info

### 2. Eigenschaften des Polyan Plus® Economic Line Blanks

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### 2. Characteristics/properties of Polyan Plus® Economic Line Blanks

- abrasion-resistant
- high breaking strength
- high flexural strength
- plaque-resistant and biocompatible

### 3. Composition

Acrylic polymer based on methyl methacrylate (PMMA)

### 4. Material-mechanical specifications

E- Modul	Ø 3.370 MPa	DIN EN ISO 20795-1:2013 compliant
Tensile strength	Ø 74,6 MPa	
Fracture stress	Ø 70,3 MPa	
Flexural strength	Ø 106,4 MPa	

### 5. Delivery forms of the Blanks

[See current CAD/CAM catalog](#)

### 6. Indication

The material is suitable for producing

- Splints

### 7. Notes on Design

Minimum wall thickness:

- 0,6 mm (zervikal)
- 1,2 mm (okklusal)

### 8. Notes on milling and further processing

We recommend using a PMMA strategy in conjunction with single-cutting or wet milling. A fine-toothed carbide bur must be used to remove the milled frameworks from the milling blank. The restorations can be finished with suitable silicone polishers and polishing brushes.

High-gloss polishing is carried out with conventional polishing pastes.

### 9. Individualization

Compounding with cold-curing polymers and composites is possible. We recommend the use of our own products or cold-curing polymers from AcrylX

### 10. Cleaning

For gentle cleaning, we recommend conventional soapy water or detergent in an ultrasonic device or hand brush.

### 11. Attachment

For the attachment of restorations made out of **Polyan Plus® Economic Line** all proven provisory and permanent cements are suitable.

### 12. Interactions

Interactions with other products commonly used in dentistry are not to be expected. If interactions cannot be ruled out in the mouth, the dental prosthesis must be removed

### 13. Storage

**Dental Plus** Blanks should be stored in a dry place at room temperature and protected from direct sunlight

### 14. Expiry date

See label on the packaging

### 15. First aid measures

Due to the properties of the material, no damage to health is normally to be expected when processing **Dental Plus** Blanks. However, should any reaction occur due to special circumstances, please observe the following instructions:

- After inhalation: Seek medical advice if symptoms occur.
- After skin contact: Rinse with water and seek medical advice if symptoms occur.
- After eye contact: In case of mechanical irritation of the eyes, rinse with water and seek medical advice.
- After ingestion: If symptoms occur, seek medical advice.

### 16. Warnings

- **The use of denture finishers can lead to embrittlement (fracture) of the material if the corresponding ingredients are used.**
- Prolonged contact with high-proof alcohol solutions should be avoided. We recommend Bacoban for disinfecting plastics.
- Please note that improper use of milling tools and polishing brushes can cause the material to overheat. This can lead to damage to the material.

0481