## Gebrauchsanweisung Dental Blank

Material:

**Dentalos Plus® mono** | **Dentalos Plus® multicolor**

### 1. Kontaktdaten

 DentalPlus GmbH  
Kohlgrub 5  
83122 Samerberg  
Tel.: 08032 9892007  
Fax: 08032 9882790  
E-Mail: info@dentalplus.info  
Internet: www.dentalplus.info

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### 2. Characteristics/properties of Dentalos Plus® Blanks

- resistant to abrasion
- high tensile strength
- high flexural strength
- plaque resistant and biocompatible
- Color fidelity and high translucency
- Residual monomer content is below 0.5% and free from benzoyl peroxide; therefore suitable for allergy sufferers
- veneering

### 3. Composition

Acrylic polymer based on methyl methacrylate

### 4. Material-mechanical specifications

E- Modul	3370 MPa	DIN EN ISO 527 compliant
Tensile strength	76,3 MPa	DIN EN ISO 527 compliant
Fracture stress	70,3 MPa	DIN EN ISO 527 compliant
Flexural strength	120 MPa	DIN EN ISO 527 compliant
Water absorption	19,36 µg/mm	DIN EN ISO 10477 compliant
Water solubility	< 1 µg/mm	DIN EN ISO 10477 compliant

### 5. Delivery forms of the Blanks

[See current CAD/CAM catalog](#)

### 6. General

**Dentalos Plus®** Blanks are milling blanks for the production of crowns and bridges for permanent use.

**Dentalos Plus®** Blanks can be used in all milling machines. Note, material overheating may be the result of improper use of cutting tools, polishing brushes, steam radiators and water baths.

This can lead to damage of the material.

### 7. Notes on Design

Minimum wall thickness:

- 0,6 mm (zervikal)
- 1,2 mm (okklusal)

### 8. Instructions for milling and processing

Milling needs to be done with a milling strategy (forward speed, number of rotation) adapted for the PMMA. Contact the manufacturer of your milling machine for further details. Do not use any diamond-coated tools, which can lead to overheating (luting) of the material.

Milling cutters, which produce good results including adequate water and air colling, are available. They are called: „**Regenbogen-Fräser**“ / „**Rainbow-Milling-Cutters**“. To uncuse the milled framework out of the milling-blank a smooth geared carbide milling cutter is used. Separate the connectors carefully and without pressure.

Restorations can be pre-polished with appropriate silicon burnishes and goat-hair-brushes. The mirror finish is carried out with Polirapid-Polishing-Paste.

To guarantee an ideal color gradient and the highest homogeneity possible, a distance of > 2 mm to the border must be adhered to.

Scrub / Ø Milling cutter	2 (2,5) mm	Scrub / Ø Milling cutter	1 mm
Rotation (rpm)	19.000 - 22.000	Rotation (rpm)	25.000 - 27.000
Feed rate x/y (mm/min)	1.800 – 2.100	Feed rate x/y (mm/min)	1.200
Feed rate z (mm/min)	1.000	Feed rate z (mm/min)	1.000 - 1.200
Z-lining	0,8 - 1,5 mm	Z-lining	0,05 - 0,08 mm
Scrub / Ø Milling cutter	2 (2,5) mm	Details	0,6 mm
Rotation (rpm)	25.000 - 27.000	Rotation (rpm)	25.000
Feed rate x/y (mm/min)	1.800 - 2.000	Feed rate x/y (mm/min)	2.000 - 2.200
Feed rate z (mm/min)	1.000	Feed rate z (mm/min)	700 - 900
Z-lining	0,2 mm	Z-lining	0,01 - 0,05
Distance to the blank outer edge	2mm	Mandatory milling geometry	1-fluted

### 9. Customizing

Custom veneers are made possible by the insertion of composites. To customize the corresponding areas are reduced so that a smooth transition to the base material is possible.

As abrasive exclusively fine carbide milling should be used. The sanded surface must be thoroughly cleaned and wetted with a bonder.

The processing of the veneering is carried out according to the manufacturer's instructions.

### 10. Cleaning

High concentrated alcohol and benzene generally lead to cracking (Crackles) in PMMA. Cleaning and disinfection fluids containing these ingredients should not be used.

Steam results in heat and pressure load, which should be avoided. Alternative cleaning possibilities are mechanical methods, or the ultrasound machine.

### 11. Attachment

For the attachment of restorations made out of **Dentalos Plus®** all proven provisory and permanent cements are suitable.

### 12. Reactions

Reactions with other products used in dentistry are not to be expected. If reactions in the mouth cannot be excluded, the denture must be removed.

### 13. Storage

Dentalos Plus® Blanks should be stored at room temperature and protected against direct sunlight.

### 14. Expiration date

check tag on the packaging.

### 15. First-Aid Measures

Due to the processing of **Dentalos Plus®** composite Blanks and its properties no health hazards are to be expected. If due to special circumstances any reactions take place, please note the following instructions:

- After inhalation: In case of problems seek medical advice.
- After skin contact: Rinse with water and seek medical advice.
- After eye contact: If mechanical irritation occurs flush eyes with water and seek medical advice.
- After ingestion: In case of discomfort seek medical advice.